



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cleanroom chair PU foam blue, height
adj. 52-65 cm, castors + foot ring**

**€434.40
plus VAT &
Shipping**

Product Images



Description

Comfortable chairs specially designed for the cleanroom, durable and functional.

Seat and back in soft blue comfort polyurethane foam (4 cm thickness), easy to decontaminate. The soft polyurethane foam is resistant to cuts, flying sparks, light alkalis, acids and lubricants.

- Seat height adjustment via gas spring
- depth adjustable
- chromed five-legged base with castors and glides
- chrome-plated pedestal, gas spring height adjustment
- Depth suspension

GS tested

Additional Information

No.	1-1166
Manufacturer (Brand)	neoLab
EAN	4058072032142
Color	silver
Material	Steel chrome plated
Chairs, stools and standing aids features	suitable for cleanroom with footrest/ring
Figure	3
Frame color	silver
Color seat / backrest	blue
Height adjustment	52 to 65 cm
Material frame	Steel chrome plated
Material rollers	Elastoplastic (TPE)
Material seat / backrest	Polyurethane (PU)
Type foot	with rollers
Type chair	Swivel / office chair

