



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## neolab® cleanroom chair PU foam blue, height adj. 55-79 cm, glider + foot ring

**Price on  
request  
plus VAT &  
Shipping**



## Description

---

### Comfortable chairs specially designed for the cleanroom, durable and functional.

Seat and back in soft blue comfort polyurethane foam (4 cm thickness), easy to decontaminate. The soft polyurethane foam is resistant to cuts, flying sparks, light alkalis, acids and lubricants.

- Seat height adjustment via gas spring
- depth adjustable
- chromed five-legged base with castors and glides
- chrome-plated pedestal, gas spring height adjustment
- Depth suspension

GS tested

## Additional Information

---

No.	1-1167
Manufacturer (Brand)	neoLab
EAN	4058072032159
Color	silver
Material	Steel chrome plated
Chairs, stools and standing aids features	suitable for cleanroom with footrest/ring
Figure	4
Height adjustment	55 to 79 cm
Material frame	Steel chrome plated
Type foot	with gliders
Frame color	silver
Color seat / backrest	blue
Material seat / backrest	Polyurethane (PU)
Type chair	Swivel / office chair

