



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® cleanroom chair PU foam black,  
height adj. 47-64 cm, with glides**

**€382.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

Durable and functional chairs especially for the cleanroom.

Seat and back are made of black PU foam with climate nubs and are therefore easy to decontaminate.

The seat height is adjustable via a gas spring, the back can be adjusted mechanically in height and inclination via fixing handles.

Five-armed, black aluminum cross base with glides.

## Additional Information

---

No.	1-1065
Manufacturer (Brand)	neoLab
EAN	4058072031725
Color	Black
Material	Aluminum anodized
Chairs, stools and standing aids features	suitable for cleanroom with adjustable seat inclination
Height adjustment	47 to 64 cm
Material frame	Aluminum anodized
Type foot	with gliders
Frame color	black
Color seat / backrest	black
Material seat / backrest	Polyurethane (PU)
Type chair	Swivel / office chair

