



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cleanroom chair PU foam black,
height adj. 47-64 cm, castors + foot ring**

**€479.70
plus VAT &
Shipping**

Product Images



Description

Durable and functional chairs especially for the cleanroom.

Seat and back are made of black PU foam with climate nubs and thus easy to decontaminate.

The seat height is adjustable via a gas spring, the back can be adjusted mechanically in height and inclination via fixing handles.

Five-armed, black aluminum cross base with soft castors and foot ring.

Additional Information

No.	1-1066
Manufacturer (Brand)	neoLab
EAN	4058072031732
Color	Black
Material	Aluminum anodized
Chairs, stools and standing aids features	suitable for cleanroom with adjustable seat inclination with footrest/ring
Height adjustment	47 to 64 cm
Material frame	Aluminum anodized
Type foot	with rollers
Frame color	black
Color seat / backrest	black
Material rollers	Elastoplastic (TPE)
Material seat / backrest	Polyurethane (PU)
Type chair	Swivel / office chair

