



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Standing aid cleanroom compatible,  
650-850 mm, disc base, seat gray**

**€444.80  
plus VAT &  
Shipping**



## Description

---

The "Labster" laboratory chair system is perfectly tailored to the needs of the laboratory in terms of form and function. It is based on the results of the Fraunhofer laboratory user study Lab I 2020, which have been consistently implemented down to the last detail.

The "Labster" standing aid is the solution for workplaces where work can only be done standing up. It enables dynamic "standing-sitting" and relieves the user's spine. With its pendulum base, it follows the work movements and opens up wide reach areas without having to stand up.

After work, "Labster" can be cleaned extremely quickly and thoroughly. Thanks to its "Hygienic Design," it has no crevices in which microorganisms, bacteria or minute particles could accumulate. All component transitions are sealed. Even the mechanics are located under a washable cover. All surfaces are resistant to all common cleaning agents and disinfectants.

The chairs are practically emission-free and meet the highest air purity classes. They are even suitable for use in cleanroom conditions.

The upholstery is washable, resistant to disinfectants and easy to clean.

The details in brief

- Air purity class 3 according to DIN EN ISO 14644-1 and class 1 according to US Fed St. 209 E
- GMP-compliant
- suitable for bio laboratories of safety classes S1, S2 and S3
- GS tested
- complies with DIN 68877
- Fraunhofer Tested® Device
- Plastic parts platinum gray, steel parts chrome-plated
- Seat height and backrest height adjustment

# Additional Information

No.	7-2014
Manufacturer (Brand)	neoLab
EAN	4058072038229
Transport temperature	Room temperature
Color	gray silver
Material	Steel chrome plated
Chairs, stools and standing aids features	cleanroom suitable ergonomically shaped
Frame color	gray silver
Color seat / backrest	gray
Height adjustment	65 to 85 cm
Material frame	Steel chrome plated
Material seat / backrest	Polyurethane (PU)
Type foot	with plate / platter
Type chair	Standing aid

