



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## neolab® Air disinfection cleaner Desi Pro 45 with UV-C radiation, for 45 sqm room area

**€1,294.00**  
**plus VAT &**  
**Shipping**



## Description

---

With the DESI **air disinfectant cleaner** you sterilize the air where and when you need it! Full offices, waiting rooms in the doctor's office, barber stores or kindergartens – the DESI air disinfectant cleaner ensures that viruses and bacteria are killed.

Depending on the ceiling height of your premises, the disinfectant cleaner **completely** cleans the air of the entire room within **5 to 8 hours!**

The DESI air disinfection cleaners are characterized by

- **NO ozone formation**
- Long life of the Philips UV lamp of up to **30 000 operating hours**
- **Complete cleaning** within a few hours
- **Quiet operation**

Especially in times of corona, it is important to protect yourself and your environment from infections. The DESI air disinfection cleaner helps you to do this! It disinfects the air via UV-C radiation from bacteria, viruses and fungal spores and differs significantly from HEPA filters. Because air disinfection with UV-C radiation:

- is more effective against microorganisms
- is more durable
- causes less follow-up costs

## Additional Information

---

No.	C-9002
Manufacturer (Brand)	neoLab
EAN	4058072160616
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	250 mm
Height	750 mm
Depth	250 mm
Power	45 W

