



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® column fan, tabletop unit, white

**€96.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **column fan** offers the optimal combination of perfect airflow and elegant design. It is equipped with 3D airflow technology, which guarantees even distribution and comprehensive oscillation. Optimal control, easy transport and numerous practical features.

Timer up to 8 h, automatic shut-off after 12 h, 3 speed settings.

Portable, with cable rewind in the base.

## Additional Information

---

No.	2-1169
Manufacturer (Brand)	neoLab
EAN	4058072158460
Transport temperature	Room temperature
Color	White
Material	Acrylonitrile butadiene styrene (ABS)
Form	around
Height	330 mm
Operation via	Membrane keypad
Number of blower stages	3
Movement	oscillating/ linear
Energy supply	Power grid
Power	30 W
Type fan	Desktop device

