



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Stakpure stainless steel ion exchanger, wall bracket

€176.00
plus VAT &
Shipping

Product Images



Description

The reliable supply of water.

Pressure-resistant ion exchange systems made of stainless steel, pressure-resistant up to 10 bar.

The economical way to fully demineralized water.

The water flows over a mixed bed of special anion and cation resins and is thus fully desalinated or demineralized by the proven method of the ion exchanger process.

Optimally applicable for the supply of ultrapure water systems, automatic analyzers, dishwashers, autoclaves or industrial use in the fields of process, battery and cooling water.

Ion exchange systems meet all current standards such as ASTM, CLSI, DIN EN 285, VDE 0510 and achieve a consistently high water quality.

Regenerations possible in our own service station!

Technical data:

- Water quality: 0,1- 20 µS/cm
- Housing: stainless steel
- Diameter: 240 mm

GSG = total salt content

Additional Information

No.	PE-0077
Manufacturer (Brand)	Stakpure
old neoLab article no.	2-0146
EAN	4066292550490
Transport temperature	Room temperature

