




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel platform trolley 820 x 545 mm, 4 swivel castors

€791.40
plus VAT &
Shipping

Product Images



Description

Universal load transporter made of stainless steel.

Hygiene-conscious design without dirt corners, with smooth sheet metal platform, beveled on all sides by 80 mm. Easy to maneuver due to 4 large swivel castors, 2 of which are lockable.

Platform: 820 x 545 mm

External dimensions:

1015 x 600 x 1020 mm (W x D x H)

Load capacity: 120 kg

Weight: 20 kg

Additional Information

No.	2-7076
Manufacturer (Brand)	neoLab
EAN	4058072008680
Transport temperature	Room temperature
Color	silver
Wide	605 mm
Height	1020 mm
Length	983 mm
Depth	983 mm
Weight	17.1 kg
Material frame	Stainless steel
Width plate	545 mm
loadable up to	120 kg
DM Rollers	125 mm
Material rollers	Elastoplastic (TPE)
Type trolley	Platform trolley

