



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® wire basket trolley with 2 baskets

**€515.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Sturdy transport cart made of plastic-coated steel tube, with two baskets.

4 casters with 10 cm Ø.

Volume:

upper basket 65 l

lower basket 24 l

Basket dimensions: 60 x 40 cm (L x W)

Load capacity per basket : 30 kg

Dimensions:

73 x 41 x 85 cm (L x W x H without handle)

Weight : 13.5 kg

## Additional Information

---

No.	2-7083
Manufacturer (Brand)	neoLab
EAN	4058072008710
Transport temperature	Room temperature
Color	Blue White
Material	Steel epoxy resin
Wide	410 mm
Width/inside	600 mm
Height	850 mm
Length	730 mm
Depth	730 mm
Weight	13.5 kg
Length inside	400 mm
Trolley features	with baskets
Number of floors	2
Number of baskets	2
DM Rollers	100 mm
Delivery	unmounted
Material frame	Steel epoxy resin
Material rollers	Elastoplastic (TPE)
Type trolley	Transport trolley
loadable up to	60 kg

