




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel transport trolley 2 shelves, 94 x 50 x 92 cm

€541.10
plus VAT &
Shipping

Product Images



Description

Trolley for transporting and depositing jars and chemicals.

As a temporary extension of the worktable for setting up experimental arrangements.

Easy to move, runs on 4 swivel castors, 2 castors with locking brake.

Hygienic, easy to clean.

Supplied fully assembled on a pallet.

Table top dimensions: 80 x 50 cm

External dimensions: 94 x 50 x 92 cm (L x W x H)

Distance between shelves: 585 mm

Load capacity: total approx. 120 kg

Weight: 15.5 kg

Additional Information

No.	2-7075
Manufacturer (Brand)	neoLab
EAN	4058072008673
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	570 mm
Height	950 mm
Length	870 mm
Depth	870 mm
Weight	15 kg
Trolley features	with special chemical resistance with parking brake
Number of floors	2
Material frame	Stainless steel
loadable up to	120 kg
DM Rollers	125 mm
Delivery	mounted
Material rollers	Elastoplastic (TPE)
Type trolley	Transport trolley

