




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel transport trolley, 4 shelves, 90 x 60 x 126 cm

€1,155.00
plus VAT &
Shipping

Product Images



Description

Transport trolley made of CrNi steel 18/9 in hygiene-conscious design.
Well maneuverable due to 4 large swivel castors, 2 castors with locking brake.
All shelves have a surrounding bead.

Shelf dimensions: 80 x 50 cm

Shelf spacing: 275 mm

Weight: 23 kg

Load capacity: 130 kg

Additional Information

No.	2-7073
Manufacturer (Brand)	neoLab
EAN	4058072008659
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	570 mm
Height	1260 mm
Length	870 mm
Depth	870 mm
Weight	22 kg
Trolley features	with special chemical resistance with parking brake
Number of floors	4
DM Rollers	125 mm
Delivery	mounted
Material frame	Stainless steel
Material rollers	Elastoplastic (TPE)
Type trolley	Transport trolley
loadable up to	160 kg

