



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® gas cylinder trolley

€335.00
plus VAT &
Shipping



Description

This **gas cylinder trolley** makes it easy for you to transport a steel gas cylinder to the desired place of use or storage. Sturdy and cleanly welded tubular steel construction with high-quality powder coating: This makes the cart permanently surface-protected, impact and scratch resistant.

Firm stand of the gas cylinder on the sheet steel bucket with three-sided border and secure hold by the hooked-in metal chain. Slipping away and tipping over not possible.

Safety handle with knuckle protection over the tubular bar wheels. With solid rubber tires, low driving noise, puncture-proof.

Suitable for bottles with 20/40 or 50 liters capacity or with a diameter of 20.4 to 22.9 cm

- Dimensions L x W: 120 x 51 cm
- Wheel diameter: 25 x 6 cm wide
- Weight: 12,5 kg
- Color: brilliant blue RAL 5007
- Load capacity 100 kg

Additional Information

No.	1-6026
Manufacturer (Brand)	neoLab
EAN	4058072036195
Color	blue
Wide	510 mm
Height	1200 mm
suitable for	Gas cylinders
Weight	12.5 kg
Number of rolls	2
loadable up to	100 kg
for bottles with volume	20 l 40 l 50 l
DM Rollers	250 mm
Material rollers	Rubber

