




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® shelf trolley supermini, 340 x 450 x 730 mm (WxDxH)

€411.30
plus VAT &
Shipping

Product Images



Description

Extra narrow shelf trolleys in top quality, ideal for small equipment or narrow spaces.

Square tubular steel bracket 20 x 20 mm, white plastic-coated.

3 panels each, lightfast plastic-coated on all sides, partially edged.

4 castors 50 mm Ø.

Additionally with 3 stabilizing rods.

Additional Information

No.	2-7041
Manufacturer (Brand)	neoLab
EAN	4058072008598
Transport temperature	Room temperature
Color	White
Material	Wood Melamine
Wide	340 mm
Height	705 mm
Length	450 mm
Depth	450 mm
Number of floors	3
Material frame	Steel epoxy resin
loadable up to	50 kg
DM Rollers	50 mm
Delivery	unmounted
Material rollers	Polycarbonate (PC)
Type trolley	Transport trolley

