



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Trolley 3 shelves, 46 x 61 x 86 cm, clear height 30 cm

€432.30
plus VAT &
Shipping

Product Images



Description

Trolley made of high-density polyethylene, ideal for laboratory and practice.

- highly loadable up to 70 kg
- electrically non-conductive
- shockproof
- easy to clean
- without sharp edges and corners
- easy assembly (rubber mallet incl.)
- Top shelf with integrated handles

4 swivel castors with 10 cm Ø (2 lockable)

basic dimensions: 46 x 61 cm (D x W)

color: gray

Additional Information

No.	2-2150
Manufacturer (Brand)	neoLab
EAN	4058072059897
Transport temperature	Room temperature
Color	gray
Material	High-density polyethylene (HDPE)
Wide	460 mm
Height	860 mm
Length	610 mm
Depth	610 mm
Distance floors	30 cm
Number of floors	2
DM Rollers	100 mm
Delivery	unmounted
Material frame	High-density polyethylene (HDPE)
Material rollers	High-density polyethylene (HDPE)
Type trolley	Transport trolley
loadable up to	90 kg

