



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Trolley 3 shelves, 46 x 61 x 86 cm, clear height 16/45 cm

**€461.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Trolley made of high-density polyethylene, ideal for laboratory and practice.

- highly loadable up to 70 kg
- electrically non-conductive
- shockproof
- easy to clean
- without sharp edges and corners
- easy assembly (rubber mallet incl.)
- Top shelf with integrated handles

4 swivel castors with 10 cm Ø (2 lockable)

basic dimensions: 46 x 61 cm (D x W)

color: gray

## Additional Information

No.	2-2151
Manufacturer (Brand)	neoLab
EAN	4058072059903
Transport temperature	Room temperature
Color	gray
Material	High-density polyethylene (HDPE)
Wide	460 mm
Height	860 mm
Length	610 mm
Depth	610 mm
Distance floors	16 cm top 45 cm bottom
Number of floors	3
DM Rollers	100 mm
Delivery	unmounted
Material frame	High-density polyethylene (HDPE)
Material rollers	High-density polyethylene (HDPE)
Type trolley	Transport trolley
loadable up to	90 kg

