




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Variocar basic trolley with 3 drawers, 60 x 45 x 73 cm

€733.70
plus VAT &
Shipping

Product Images



Description

Mobile side table with worktop, drawers and storage compartment at the bottom.
Drawers with smooth running ball telescopic runners with locking mechanism.

- chrome-plated square tubular steel frame
- Tops and drawers with high-quality, five-layer coating on both sides
- 4 smooth-running castors with 75 mm Ø
- Total load up to approx. 50 kg

The trolley is delivered disassembled, but can be easily assembled.

Dimensions: 600 x 450 x 730 mm (W x D x H)

Other models on request.

Additional Information

No.	2-7022
Manufacturer (Brand)	neoLab
EAN	4058072008482
Transport temperature	Room temperature
Color	silver white
Material	Wood Melamine
Wide	450 mm
Height	730 mm
Length	600 mm
Depth	600 mm
Trolley features	with drawers
Number of drawers	3
Material frame	Steel chrome plated
loadable up to	50 kg
DM Rollers	75 mm
Delivery	unmounted
Material rollers	Polycarbonate (PC)
Type trolley	Transport trolley Support/storage table

