



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® transport trolley, 2 shelves

€235.50
plus VAT &
Shipping

Product Images



Description

For transport and parking. 4 swivel castors, 2 with locking brake.

Made of stainless steel, hygienic, easy to clean.

Space-saving delivery for self-assembly.

Shelf dimensions: 800 x 500 mm

External dimensions: 950 x 500 x 950 mm (LxWxH)

Load capacity per shelf: max. 40 kg

Additional Information

No.	E-7075
Manufacturer (Brand)	neoLab
EAN	4058072015008
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	500 mm
Width/inside	500 mm
Height	950 mm
Length	950 mm
Depth	950 mm
Trolley features	with parking brakeelectrically conductive
Distance floors	58.5 cm
Number of floors	2
DM Rollers	125 mm
Weight	15.5 kg
Length inside	800 mm
Material frame	Stainless steel
Material rollers	Polyamide (PA. non-specific)
Type trolley	Transport trolley
loadable up to	100 kg

