




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Variocar undercarriage, 3 drawers, 60 x 45 x 52 cm

€761.40
plus VAT &
Shipping

Product Images



Description

Low trolley with an overall height of only 52 cm.

Fits under the desk or examination couch.

- chrome-plated steel frame
- high quality melamine resin coating of the tops and drawers
- lightfast, abrasion resistant, resistant to many chemicals
- Drawers with smooth-running telescopic ball bearing runners and locking mechanism
- 4 castors with 50 mm Ø

The trolley is delivered disassembled, but is easy to assemble.

Dimensions: 60 x 45 x 52 cm (L x W x H)

Additional Information

No.	2-7031
Manufacturer (Brand)	neoLab
EAN	4058072008543
Transport temperature	Room temperature
Color	silver white
Material	Wood Melamine
Wide	450 mm
Height	520 mm
Length	600 mm
Depth	600 mm
Trolley features	with special chemical resistance with drawers
Number of drawers	3
Material frame	Steel chrome plated
loadable up to	50 kg
DM Rollers	50 mm
Delivery	unmounted
Material rollers	Polycarbonate (PC)
Type trolley	Transport trolley Support/storage table

