



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Rotor Blender basic unit without attachment

**€1,723.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Blender-homogenizer with built-in timer from 1 to 60 seconds and indicator light.

Robust device with powerful electric motor with ball bearings.

The stainless steel blades are dynamically balanced and equipped with self-lubricating, maintenance-free bearings.

The base is enameled.

With 2-step switch for 18,500 and 21,000 rpm with pilot lamp.

Operation with various attachments.

Weight basic unit: 4 kg

Dimensions: 22 cm Ø, 52 cm height

Power supply:

220 V, 50 Hz

## Additional Information

---

No.	3-2095
Manufacturer (Brand)	neoLab
EAN	4058072018917
Transport temperature	Room temperature
Speed MAX	22000 rpm
Color	White
Height	520 mm
Blender and mixer features	with timer with continuous function mixing attachment available as accessory
Weight	4 kg
Voltage	230 V
Number of speed levels	2
Functionality	electric
Power	860 W
Type homogenizer	Stand homogenizer/blender

