




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## WLD-TEC Annealing carousel for inoculation loops Autoloop PRO

**€1,229.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Fully automatic eyelet carousel<sup>autoloopPRO</sup> made of stainless steel, including 4 inoculation loop holders. By simply plugging it on, any Fuego SCS burner becomes an annealing station.

The fully graphic display shows flame time, start options, cooling time, position of cooled inoculation loops, sterilization monitor.

Fully automatic control of the annealing process, quick start button.

Stainless steel, anodized aluminum housing.

Automatic switch-off after 1 hour.

Overtemperature protection, holder removal monitoring, flame activation indicator.

Operation via push button or foot pedal.

Glow time: 1-15 sec. adjustable

Dimensions base plate: 135 x 180 (W x D)

**Height with carousel:** 305 mm

Weight: 1200 g

Other models for other burners and a wide range of accessories on request.

## Additional Information

---

No.	WT-0097
Manufacturer (Brand)	WLD-TEC
old neolab article no.	7-8205
EAN	4058072443993
Transport temperature	Room temperature
Material	Stainless steel
Wide	180 mm
Height	305 mm
Length	135 mm
Weight	1.2 kg
Number of holders	4

