



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC Sterimax electric annealing system

€601.00
plus VAT &
Shipping

Product Images



Description

Annealing **without flame !**

The Sterimax electric annealing system is ideally suited for annealing inoculation loops, needles and instruments. Ideal for laboratories and safety workbenches where gas and open flames must not be used.

The annealing tube made of wear-resistant quartz glass reaches the optimum annealing temperature of 900°C (1650°F) after only a few minutes.

For annealing, simply hold the inoculating loop in the quartz glass tube.

The inoculating loop is sterilized after only 5-7 seconds.

With integrated splash guard. The stainless steel housing is also suitable for extreme laboratory conditions. The unit has a 7-position angle adjustment and has a low, stable housing. A residual heat indicator is integrated as a safety feature.

Additional Information

No.	WT-0033
Manufacturer (Brand)	WLD-TEC
old neoLab article no.	7-8215
EAN	4058072443900
Transport temperature	Room temperature
Temp MAX	950 °C
Temp MIN	900 °C
Wide	175 mm
Height	129 mm
Length	139 mm
Depth	139 mm
Type of sterility	heat sterile
Operation via	Knobs/Buttons
Functionality	electric
Weight	0.96 kg
Power	160 W

