



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC Sterimax smart, infrared sterilizer for inoculation loops

€959.00
plus VAT &
Shipping

Product Images



Description

The Sterimax smart sterilizes your inoculating loops and similar small instruments quickly and efficiently using IR radiation. No heating phase or gas cartridges are required. The device is ready for immediate use. Your instruments are sterilized in a special quartz glass annealing tube in 5 to 10 seconds.

- Temperatures: 750 to 1000 °C (outer housing is not heated)
- Annealing tube can be easily removed for cleaning
- Operation via sensitive IR sensor touch function
- with individually adjustable timer function for 2 users
- with visual and acoustic signal at the end of the process

Additional Information

No.	WT-0021
Manufacturer (Brand)	WLD-TEC
old neoLab article no.	7-8214
EAN	4058072443894
Transport temperature	Room temperature
Temp MAX	1000 °C
Temp MIN	750 °C
Wide	170 mm
Height	180 mm
Length	110 mm
Depth	110 mm
Operation via	Membrane keypad
Weight	1.8 kg
Type of sterility	radiation sterile
Functionality	electric electronic
Power	600 W

