



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC laboratory gas burner Powerjet

€272.00
plus VAT &
Shipping

Product Images



Description

The solution for all flame protection sterilizations in microbiology or biotechnology.

One press of the button and the burner is ignited. The flame can be continuously adjusted in size and goes out when the button is released.

In addition, both Powerjets are equipped with a locking device for fixing the flame. This allows mobile short-time or continuous flaming as required.

Order No. 7-8227

The Powerjet is suitable for CGI750 gas cartridges (liquid gas).

Flame temperature: 1350°C

Nominal heat load: 1.54 kW

Dimensions : 150 x 180 mm (W x H)

Weight: 300 g

Additional Information

No.	WT-0089
Manufacturer (Brand)	WLD-TEC
old neoLab article no.	7-8227
EAN	4058072443979
Transport temperature	Room temperature
Wide	150 mm
Height	180 mm
Display (reading, setting)	analog
Operation via	Knobs/Buttons
Weight	0.3 kg
Energy supply	without power supply
Gas type	Liquid gas (propane. Butane)
Power	1540 W
Gas supply via	Gas cartridge
Model burner	Powerjet (7-8227)
Type burner	Cartridge burner

