



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC Laboratory gas burner Powerjet2 w. Gas cartridge adapter & stainless steel holder

€448.00
plus VAT &
Shipping

Product Images



Description

The solution for all flame protection sterilizations in microbiology or biotechnology.

One press of the button and the burner is ignited. The flame can be continuously adjusted in size and goes out when the button is released.

In addition, both Powerjets are equipped with a locking device for fixing the flame. This allows mobile short-time or continuous flaming as required.

Order No. 7-8229

ThePowerjet2 is suitable for natural gas and propane/butane gas.

Included in delivery are a gas cartridge adapter for CV360 cartridges, a stainless steel holder as well as nozzles for natural gas and propane/butane gas.

Flame temperature: 1350°C for LPG, 1300°C for natural gas

Rated heat load: 1.7 kW for LPG, 1 kW for natural gas

Dimensions: 150 x 210 mm (W x H)

Weight: 500 g

Additional Information

No.	WT-0090
Manufacturer (Brand)	WLD-TEC
old neolab article no.	7-8229
EAN	4058072443986
Transport temperature	Room temperature
Wide	150 mm
Height	210 mm
Burner properties	with continuous function
Display (reading, setting)	analog
Operation via	Knobs/Buttons
Weight	0.5 kg
Energy supply	without power supply
Gas type	Natural gas Liquid gas (propane. Butane)
Power	1700 W
Gas supply via	Gas cartridge
Model burner	Powerjet 2 (7-8229)
Type burner	Cartridge burner

