




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hygiene PC mouse, large, 5 buttons, IP 69

€167.80
plus VAT &
Shipping

Product Images



Description

The **neolab Hygiene Silicone PC Mouse** is completely sealed so it can be wiped or sprayed and even dipped in a standard surface disinfectant. It is the ideal pointing device for operating rooms, clean rooms, laboratories or intensive care units, as well as anywhere fluid incidents need to be considered.

Manufactured with advanced technology, the Hygiene PC Mouse comes with EN60601 and **IP68 certifications**. The Hygiene PC Mouse is an optical mouse with 5-button design for scrolling and clicking. This full size mouse provides a good feel for best user experience.

Additional Information

No.	1-0120
Manufacturer (Brand)	neolab
EAN	4058072188641
autoclavable at 121°C	No
Color	White
sterile	No
Weight	142 g
Interface	USB
Protection class	IP 68
for operating system	Windows Linux Mac OS
Type computer accessories	USB hygiene mouse large

