



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA Digital Imaging and Analysis System L-340

**Price on
request
plus VAT &
Shipping**

Product Images

Description

SERVA's Digital Imaging and Analysis System III is a fast and convenient gel documentation system that masters the daily tasks of documentation and 1D gel analysis. Solid hardware combined with a digital SLR camera and easy to understand Labimage L-340 1D analysis software provides an excellent tool to meet your needs. The instrument is equipped with a UV filter and filter holder. UV, blue and white light transilluminator or epi-white light are optional. The easy to understand L340 1D analysis software allows for advanced 1D analysis.

Components:

- Darkroom (approx. 42 x 55 x 52 cm)
- UV filter and filter holder
- Digital SLR camera
- Labimage L-340

Additional Information

No.	SV-0205
Manufacturer (Brand)	SERVA
EAN	4058072373429
Transport temperature	Room temperature
Wide	5500 mm
Height	5200 mm
Depth	4200 mm
Weight	12 kg
EPI light source	optional UV/BLEND/WL
Filter	580 nm

