



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Electrophoresis Chamber Maxi horizontal for gel size 204 x 250 mm + 204 x 200 mm

**€647.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Electrophoresis chamber Maxi horizontal for gel size 204 x 250 mm and 204 x 200 mm.

Including gel carrier for gel size 204 x 250 mm and gel casting stand

Buffer volume 1150 ml

Dimensions: 37 x 26.5 x 9.5 cm (L x W x H)

## Additional Information

---

No.	7-0219
Manufacturer (Brand)	neolab
EAN	4058072027759
Transport temperature	Room temperature
Length	37 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand gel carrier is UV-permeable
Number of comb positions	4
Alignment	horizontal
Wide chamber	26.5 cm
Size chamber	Maxi
Height chamber	9.5 cm
Dimensions gel	200 x 204 mm 204 x 250 mm
Separation lines	4 x 6 cm 3 x 8 cm 2 x 12 cm 1 x 24 cm
Volume gel	225 ml
Volume chamber	1150 ml

