



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Electrophoresis Chamber Maxi-wide horizontal for gel size 239 x 145 mm

€644.70
plus VAT &
Shipping

Product Images



Description

Maxi-wide horizontal for gel size 204 x 250 mm and 204 x 200 mm

incl. gel carrier for gel size 204 x 250 mm and gel casting stand

Buffer volume 1050 ml

Dimensions: 29.5 x 29.5 x 9.5 cm (L x W x H)

Additional Information

No.	7-0234
Manufacturer (Brand)	neoLab
EAN	4058072027896
Transport temperature	Room temperature
Length	29.5 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand gel carrier is UV-permeable
Number of comb positions	4
Height chamber	9.5 cm
Dimensions gel	145 x 239 mm
Alignment	horizontal
Wide chamber	29.5 cm
Size chamber	Maxi
Separation lines	4 x 6 cm 3 x 8 cm 2 x 12 cm 1 x 24 cm
Volume gel	225 ml

