



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Midi-I Bufferflow

**€554.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Midi-I Bufferflow

incl. gel carrier for gel size 151 x 170 mm and gel casting stand.

Buffer volume 850 ml

Dimensions: 29.5 x 21.0 x 9.5 cm (L x W x H)

## Additional Information

---

No.	7-0215
Manufacturer (Brand)	neoLab
EAN	4058072027742
Transport temperature	Room temperature
Length	29.5 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand suitable for use with magnetic stirrer gel carrier is UV-permeable
Number of comb positions	3
Height chamber	9.5 cm
Dimensions gel	100 x 151 mm 151 x 170 mm
Alignment	horizontal
Wide chamber	21 cm
Size chamber	Midi
Separation lines	3 x 5.3 cm 2 x 8.0 cm 1 x 16 cm
Volume gel	128 ml
Volume chamber	850 ml

