



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Electrophoresis Chamber Mini horizontal for gel size 81 x 100 mm + 81 x 70 mm

€430.40
plus VAT &
Shipping

Product Images



Description

Electrophoresis chamber Mini horizontal for gel size 81 x 100 mm and 81 x 70 mm.

Including gel carrier for gel size 81 x 100 mm and gel casting stand.

Buffer volume 400 ml.

Dimensions: 22 x 13.5 x 9.5 cm (L x W x H)

Additional Information

No.	7-0155
Manufacturer (Brand)	neolab
EAN	4058072027292
Transport temperature	Room temperature
Length	22 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand gel carrier is UV-permeable
Number of comb positions	2
Alignment	horizontal
Wide chamber	13.5 cm
Size chamber	Mini
Height chamber	9.5 cm
Dimensions gel	81 x 100 mm 70 x 81 mm
Separation lines	2 x 4.5 cm 1 x 9.5 cm
Volume gel	43 ml
Volume chamber	400 ml

