



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Axygen® horizontal electrophoresis device, 15 cm

€737.82
plus VAT &
Shipping

Product Images



Description

Horizontal gel chambers for agarose gel electrophoresis.

Included in delivery:

One chamber, one cover and electrodes.

Separate gel casting stand; this allows you to separate your samples and cast another gel at the same time.

Gel carriers in 2 different sizes each: 15 x 10 cm, 15 x 15 cm

Combs: 1 x 16 compatible multichannel teeth, thickness: 1.5 mm
and 1 x 20 teeth, thickness: 1.5 mm

Contrast plate for easier sample application.

The electrodes can be easily replaced.

Additional Information

No.	CQ-0413
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2002
EAN	4058072226022
Transport temperature	Room temperature
Length	15 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand suitable for use with magnetic stirrer gel carrier is UV-permeable
Number of comb positions	4
Height chamber	3.3 cm
Dimensions gel	150 x 150 mm
Alignment	horizontal
Wide chamber	15 cm
Size chamber	Midi
Volume gel	225 ml

