



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Axygen® horizontal electrophoresis instrument, 20 cm

**€827.42**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Horizontal gel chambers for agarose gel electrophoresis.

Included in delivery:

One chamber, one cover and electrodes.

Separate gel casting stand; this allows you to separate your samples and cast another gel at the same time.

Gel carriers in 2 different sizes each: 20 x 20 cm, 20 x 25 cm

**Combs:** 2 x 20 compatible multichannel teeth, thickness: 1.5 mm  
and 2 x 40 compatible multichannel teeth, thickness: 1.5 mm

Contrast plate for easier sample application.

The electrodes can be easily replaced.

## Additional Information

---

No.	CQ-0444
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2003
EAN	4058072226169
Transport temperature	Room temperature
Length	25 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand suitable for use with magnetic stirrer gel carrier is UV-permeable
Number of comb positions	4
Alignment	horizontal
Wide chamber	20 cm
Size chamber	Maxi
Height chamber	3.5 cm
Dimensions gel	200 x 250 mm
Volume gel	400 ml

