



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Axygen® horizontal electrophoresis device, 7 cm

€413.67
plus VAT &
Shipping

Product Images



Description

Horizontal gel chambers for agarose gel electrophoresis.

Included in delivery:

One chamber, one cover and electrodes.

Separate gel casting stand; this allows you to separate your samples and cast another gel at the same time.

Gel carriers in 2 different sizes each: 7 x 7 cm, 7 x 10 cm

Comb: 2 x 8 teeth, thickness: 1.5 mm

Contrast plate for easier application of samples.

The electrodes can be easily exchanged.

Additional Information

No.	CQ-0459
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2000
EAN	4058072226282
Transport temperature	Room temperature
Length	10 cm
Electrophoresis Equipment Features	with gel carrier with gel casting stand gel carrier is UV-permeable
Number of comb positions	4
Height chamber	3 cm
Dimensions gel	70 x 100 mm
Alignment	horizontal
Wide chamber	7 cm
Size chamber	Mini
Volume gel	70 ml

