



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA SERVAGel™ Neutral HSE, 2D ready-to-use gel, 2D sample bag

€158.00
plus VAT &
Shipping

Product Images



Description

The SERVAGel™ Neutral HSE ready-to-use gel is a version of SERVAGel™ Neutral pH 7.4 specifically designed for high-speed electrophoresis. Two 2D gels can be run with the standard SDS Tris glycine buffer system (according to Laemmli) at 300 V, reducing electrophoresis time to 40 min

The separation range is between 6.5 up to 200 kDa. For the first dimension, an IPG strip of 7 cm length can be used. The separation distance of the gel is 7 cm. The in-house developed neutral buffer system guarantees a significantly extended shelf life compared to classical Tris-Glycine gels with a Laemmli buffer system.

Additional Information

No.	SV-0114
Manufacturer (Brand)	SERVA
EAN	4058072372514
Transport temperature	Wet ice
Number of wells	1
Thickness	1 mm

