



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA SERVAGel™ Neutral HSE gel ready-to-use, 12 sample bags

€158.00
plus VAT &
Shipping

Product Images



Description

The SERVAGel™ Neutral HSE ready-to-use gel is a version of the SERVAGel™ Neutral pH 7.4 Gradient developed specifically for high-speed electrophoresis. The gel can be run with the standard SDS-Tris glycine buffer system (according to Laemmli) at 400 V, which reduces the electrophoresis time to 20 min.

The separation range is between 6.5 up to 200 kDa. The in-house developed neutral buffer system guarantees a significantly extended shelf life compared to classical Tris-glycine gels with a Laemmli buffer system.

Additional Information

No.	SV-0110
Manufacturer (Brand)	SERVA
EAN	4058072372477
Transport temperature	Wet ice
Number of wells	12
Thickness	1 mm

