



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SERVA SERVAGel™ Neutral HSE gel ready-to-use, 15 sample bags

Price on  
request  
plus VAT &  
Shipping

### Product Images



## Description

---

The SERVAGel™ Neutral HSE ready-to-use gel is a version of the SERVAGel™ Neutral pH 7.4 Gradient developed specifically for high-speed electrophoresis. The gel can be run with the standard SDS-Tris glycine buffer system (according to Laemmli) at 400 V, which reduces the electrophoresis time to 20 min.

The separation range is between 6.5 up to 200 kDa. The in-house developed neutral buffer system guarantees a significantly extended shelf life compared to classical Tris-glycine gels with a Laemmli buffer system.

## Additional Information

---

No.	SV-0117
Manufacturer (Brand)	SERVA
EAN	4058072372545
Transport temperature	Wet ice
Number of wells	15
Thickness	1 mm

