



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA SERVAGel™ Neutral pH 7.4 Gradient Gel Ready-to-use Gel, 15 Sample Bags

**Price on
request
plus VAT &
Shipping**

Product Images



Description

The SERVAGel™ Neutral pH 7.4 Gradient Finishing Gel has all the advantages of SERVAGel™ Neutral pH 7.4 , but has even better resolution (from 5 kDa up to approximately 100 kDa) when separating low molecular weight proteins. The separation range is from 5 kDa up to 200 kDa. To improve the resolution and separation sharpness of the protein bands, SERVA Tris-Tricin running buffer and Tris-Tricin/SDS sample buffer should be used.

Additional Information

No.	SV-0120
Manufacturer (Brand)	SERVA
EAN	4058072372576
Transport temperature	Wet ice
Number of wells	15
Thickness	1 mm

