



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SERVA SERVAGel™ TG PRiME™ 14 % ready-to-use gel, 15 sample bags

Price on  
request  
plus VAT &  
Shipping

### Product Images



## Description

---

The SERVAGel™ TG PRIME™ 14 % ready-to-use gel has an in-house developed modified Tris/glycine buffer system, which not only guarantees a significantly extended shelf life compared to classical Tris/glycine gels according to Laemmli, but also allows short electrophoresis times when using standard Tris/glycine sample and run buffers (2 gels, 35 min at 250 V const.). It can be used under native and denaturing conditions. The separation range is from 10 up to 100 kDa.

## Additional Information

---

No.	SV-0163
Manufacturer (Brand)	SERVA
EAN	4058072373009
Transport temperature	Wet ice
Number of wells	15
Thickness	1 mm

