



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Axygen® gel carrier 7 x 10 cm for use with 7 cm gel box, UV transparent

€66.39
plus VAT &
Shipping

Product Images



Description

Gel carrier for electrophoresis chamber Axygen®, horizontal **Size: 10 x 10 cm**

Additional Information

No.	CQ-0395
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2013
EAN	4058072225926
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	100 mm
Height	33 mm
Length	100 mm
Depth	100 mm
Electrophoresis Equipment Features	with gel carrier
Number of comb positions	4
Dimensions gel	100 x 100 mm
Volume gel	100 ml

