



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Axygen® gel carrier 15 x 10 cm for use with 15 cm gel box, UV transparent

€73.13
plus VAT &
Shipping

Product Images



Description

Gel carrier in for electrophoresis chamber Axygen®, horizontal **Size: 15 x 10 cm**

Additional Information

No.	CQ-0416
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2014
EAN	4058072226039
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	100 mm
Height	33 mm
Length	150 mm
Depth	150 mm
Electrophoresis Equipment Features	with gel carrier
Number of comb positions	4
Dimensions gel	100 x 150 mm
Volume gel	150 ml

