



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Axygen® gel carrier 15 x 15 cm for use with 15 cm gel box, UV transparent

**€73.13**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Gel carrier for electrophoresis chamber Axygen®, horizontal **Size: 15 x 15 cm**

## Additional Information

---

No.	CQ-0422
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2015
EAN	4058072226060
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	150 mm
Height	33 mm
Length	150 mm
Depth	150 mm
Electrophoresis Equipment Features	with gel carrier
Number of comb positions	4
Dimensions gel	150 x 150 mm
Volume gel	225 ml

