



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Axygen® gel carrier 20 x 20 cm for use with 20 cm gel box, UV transparent

**€100.45**  
plus VAT &  
Shipping

### Product Images



## Description

---

Gel carrier for electrophoresis chamber Axygen®, horizontal **Size: 20 x 20 cm**

## Additional Information

---

No.	CQ-0446
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2016
EAN	4058072226176
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	200 mm
Height	35 mm
Length	200 mm
Depth	200 mm
Electrophoresis Equipment Features	with gel carrier
Number of comb positions	4
Dimensions gel	200 x 200 mm
Volume gel	400 ml

