



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Axygen® gel carrier 7 x 10 cm for use with 7 cm gel box, UV transparent

€57.17
plus VAT &
Shipping

Product Images



Description

Gel carrier for electrophoresis chamber Axygen®, horizontal **Size: 7 x 10 cm**

Additional Information

No.	CQ-0462
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2011
EAN	4058072226312
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	70 mm
Height	30 mm
Length	100 mm
Depth	100 mm
Electrophoresis Equipment Features	with gel carrier
Number of comb positions	4
Dimensions gel	70 x 100 mm
Volume gel	50 ml

