



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SAF Wärmetechnik LabHeat heating mantle for round bottom flasks 500 ml

**€321.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

LabHeat standard heating mantles for round bottom flasks in excellent quality, hand braided.

- 50 and 100 ml with 1 heating zone
- from 250 ml with 2 heating zones (power optionally adjustable in 1/4 steps)

Heating conductor temperature: max. 450°C

Accessories: metal tripod (Order No. 2-3950 ff.)  
for all sizes, except 50 ml.

## Additional Information

---

|                                  |                    |
|----------------------------------|--------------------|
| No.                              | WF-0422            |
| Manufacturer (Brand)             | SAF Wärmetechnik   |
| old neoLab article no.           | 3-2228             |
| EAN                              | 4058072428174      |
| Transport temperature            | Room temperature   |
| Color                            | White              |
| Material                         | Glass silk         |
| suitable for                     | Round bottom flask |
| Voltage                          | 230 V              |
| Functionality                    | electric           |
| Power                            | 250 W              |
| for vessels with volume up to    | 500 ml             |
| Heatsealing band temperature MAX | 450 °C             |
| for vessels with DM up to        | 105.0 mm           |

