




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SAF Wärmetechnik LabHeat heating mantle for round bottom flasks 1000 ml

€360.00
plus VAT &
Shipping

Product Images



Description

LabHeat standard heating mantles for round bottom flasks in excellent quality, hand braided.

- 50 and 100 ml with 1 heating zone
- from 250 ml with 2 heating zones (power optionally adjustable in 1/4 steps)

Heating conductor temperature: max. 450°C

Accessories: metal tripod (Order No. 2-3950 ff.)
for all sizes, except 50 ml.

Additional Information

No.	WF-0423
Manufacturer (Brand)	SAF Wärmetechnik
old neoLab article no.	3-2229
EAN	4058072428181
Transport temperature	Room temperature
Color	White
Material	Glass silk
suitable for	Round bottom flask
Voltage	230 V
Functionality	electric
Heatsealing band temperature MAX	450 °C
Power	450 W
for vessels with DM up to	131.0 mm
for vessels with volume up to	1000 ml

