



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SAF Wärmetechnik LabHeat heating mantle for round bottom flasks 5 l

**€658.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

LabHeat standard heating mantles for round bottom flasks in excellent quality, hand braided.

- 50 and 100 ml with 1 heating zone
- from 250 ml with 2 heating zones (power optionally adjustable in 1/4 steps)

Heating conductor temperature: max. 450°C

Accessories: metal tripod (Order No. 2-3950 ff.)  
for all sizes, except 50 ml.

## Additional Information

---

No.	WF-0428
Manufacturer (Brand)	SAF Wärmetechnik
old neoLab article no.	3-2232
EAN	4058072428211
Transport temperature	Room temperature
Color	White
Material	Glass silk
suitable for	Round bottom flask
Voltage	230 V
Functionality	electric
Power	1200 W
for vessels with volume up to	5000 ml
Heatsealing band temperature MAX	450 °C
for vessels with DM up to	223.0 mm

