



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® automatic hotplate 500 W, plate 80 mm Ø

**€75.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The smallest automatic hotplate in the world.

Stainless steel plastic housing with a heating surface of 80 mm diameter. Ideal for heating small vessels.

- automatic temperature control
- Heating plate infinitely variable up to 350 °C
- Overheating protection
- Cable rewind in the base housing
- 14 x 14 x 7 cm high

Power: 500 Watt

Connection: 230 Volt/50 Hz

## Additional Information

---

No.	2-3869
Manufacturer (Brand)	neoLab
EAN	4058072001629
Transport temperature	Room temperature
Color	black silver
Material	Cast iron
Temp MAX	350 °C
Wide	140 mm
Height	70 mm
Length	140 mm
Depth	140 mm
Operation via	Knobs/Buttons
Voltage	230 V
Heaters and warming cabinets features	with stepless temperature control
Functionality	electric
Power	500 W
Material housing	Stainless steel polyamide (PA. non-specific) glass fiber reinforced

