



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® flash cooking plate, 180 mm Ø

€141.20
plus VAT &
Shipping

Product Images



Description

Inexpensive flash hotplate, also for laboratory use.
With 7-level power control. Housing white enameled.
Hotplate: 180 mm Ø

Technical data

Power: 1,500 watts

Connection: 230 V/50 Hz

Dimensions: 28.5 x 31.5 x 8 cm

Weight: 3.3 kg

Additional Information

No.	2-3871
Manufacturer (Brand)	neoLab
EAN	4058072001636
Transport temperature	Room temperature
Color	White
Material	Cast iron
Temp MAX	260 °C
Wide	315 mm
Height	80 mm
Length	285 mm
Depth	285 mm
Operation via	Knobs/Buttons
Weight	3.3 kg
Voltage	230 V
Functionality	electric
Power	1500 W
Material housing	Enamel

