



neoLab® induction hotplate, 2000 W

Product Images



Description

The fastest and most effective way to cook!

- no heat-up time, immediate power transfer
- Schott Ceran® cooking surface, easy to clean
- Heating zone 80 - 200 mm Ø
- Power selection in 9 levels
- Cooking time preselection, automatic switch-off
- Residual heat indicator

Model 2-3874 Induction hotplate

- electron. Control via touch control sensors
- Temperature preselection in 10 steps (60° to 240°C)
- 2000 W
- Dimensions: 330 x 380 x 60 mm
- Cooking zone 100 - 200 mm Ø

Additional Information

No.	2-3874
Manufacturer (Brand)	neoLab
EAN	4058072001650
Color	silber anthrazit
Material	Ceran - Glas
Temp MAX	240 °C
Wide	330 mm
Height	60 mm
Length	380 mm
Operation via	Sensor
Weight	3.2 kg
Voltage	230 V
Functionality	elektronisch elektrisch
Power	2000 W

